

chunc spica

post operative mobility



www.chunc.com



- Rental service available – visit www.chunc.com
 - Suitable for indoor and outdoor use
 - Individually removable adjustable leg splints
 - Three sizes: xs, small, medium
 - Set back recline angle between 90° and 160°
 - User mass: 54kg xs, small, medium
 - 20° posterior, 2° anterior tilt
 - Seat depth growth of 7.5cm/3" in each size
- Suitable as a forward-facing seat in a motor vehicle
- Crash tested ISO 7176 (ANSI RESNA) part 19

The Chunc Spica offers a post operative mobility solution for young people in hip spica casts and other requirements such as radial leg extension braces. Highly configurable leg supports and backrest are able to accommodate the awkward position created by the large post operative plaster.

The Chunc Spica facilitates early discharge from the hospital environment to the comfort and security of the home. Furthermore, having being crash tested (suitable for transportation) the chair allows continued participation in the social, education and therapy environments.

Support Options

Range of Headrests



Thoracic Support



Spica Wing Support



Range of Harnesses



Support Options (continued)

Multi-plane Armrests
PSIS Insert


Seating Options

<i>Basic Seat</i>	<i>Basic Backrest</i>
<i>Contour Seat</i>	<i>Supportive Backrest</i>
<i>Ischial Seat</i>	<i>Enhanced Contour Seat</i>

Other Options

<i>Running Disc Brakes</i>	<i>Spoke Guards</i>
<i>Tray</i>	<i>Rain Cover</i>
<i>Comm. Aid Mount</i>	<i>Poncho</i>

Colour Options



chunc spica

Metric

Imperial

Sizing Guide	
X-SMALL	
A	36-44 (14.1-17.3)
B	21-43 (8.2-16.9)
C	24-30 (9.4-11.9)
D	15-26 (5.9-10.2)
SMALL	
A	42-52 (16.5-20.4)
B	21-43 (8.2-16.9)
C	30-36 (11.8-14.3)
D	20-31 (7.9-12)
MEDIUM	
A	48-60 (18.8-23.6)
B	21-43 (8.2-16.9)
C	35-41 (13.8-16.3)
D	25-36 (9.8-14)

Product specification	X-Small	Small	Medium	X-Small	Small	Medium
WEIGHT						
Maximum user weight	54kg	54kg	54kg	119lb	119lb	119lb
Product weight	24kg	24.25kg	24.5kg	52.8lb	53.35lb	53.9lb
Chassis weight	16.25kg	16.5kg	16.75kg	35.75lb	36.3lb	36.85lb
DIMENSIONS						
Product height	126	126	126	49.5	49.5	49.5
Dismantled product L, W, H	82,56,54	82,56,54	82,56,54	32,22,21	32,22,21	32,22,21
Front castor diameter	18	18	18	7	7	7
Rear wheel diameter	38	38	38	15	15	15
SEATING DIMENSIONS						
Seat cushion depth range	24-30	30-36	35-41	9.4-11.9	11.8-14.3	13.8-16.3
Seat cushion depth growth	6	6	6	2.5	2.5	2.5
Seat width	26	31	36	10.2	12.2	14.2
Back rest height range	36-44	42-52	48-60	14.2-17.3	16.5-20.4	18.8-23.6
Back rest height growth	8	10	12	3.1	4	4.7
Back rest width	21	25	29	8.2	9.8	11.4
Leg splint angle range	0-20°	0-20°	0-20°	0-20°	0-20°	0-20°
Leg splint length range	21-38	21-43	21-48	8.3-15	8.3-16.9	8.3-18.9
Hip abduction range	0-45°	0-45°	0-45°	0-45°	0-45°	0-45°
Back rest recline angle range	0 to -70°	0 to -70°	0 to -70°	0 to -70°	0 to -70°	0 to -70°
Tilt in space	+2 to -20°	+2 to -20°	+2 to -20°	+2 to -20°	+2 to -20°	+2 to -20°
	<i>Dimensions in cm</i>			<i>Dimensions in inches</i>		

HRS document number 1101-03-08

